HEALTHCARE SOLUTION™

Body Fat Analyzer

GBF-1257

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Thank you for your purchase.

Please do read the user manual carefully and thoroughly so as to ensure the safe usage of this product, and keep the manual well for further reference in case you have problems.

TABLE OF CONTENT

Safety Information

Safety and Usage Information	23
Tips	24

Overview

Device Components	24
Function Indication	25
List	25
LCD Display	25

Initial Start-Up

General Instructions	26
Insert the Batteries	26
Set Up Your Profile	26

Start Measuring

Initialising Your Scale	27
Weight Only Operation	28
First Measurement	28
Daily Measurement	29
BIA Measurement	30

Manage Your Health

Body Fat Indicator	31
--------------------------	----

Troubleshooting

Error Prompt	31
When Measuring	32

Specifications	33
----------------------	----

Maintenance	34
-------------------	----

Warranty 34

Appendix

The body fat ratio 34

The body muscle ratio 34

The bone range 35

Health Tips – About Body Fat 35

EMC Guidance 35

SAFETY AND USAGE INFORMATION

The warning signs and symbols are essential to ensure your correct and safe use of this product and protect you and others from injury. Please kindly find the meanings of the warning signs and symbols, which you may encounter in the label and user manual, as follows:

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "COMPLIES WITH EU REQUIREMENTS"
	Symbol for "MANUFACTURE DATE"		Symbol for "MANUFACTURER"
SN	Symbol for "SERIAL NUMBER"		Symbol for "ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "DIRECT CURRENT"		

CAUTION

CAS Body Fat Analyzer GBF-1257 offers you a seamless way to manage your health. Please be aware that this device is designed for healthy population over 18 years old self-measuring and self-monitoring body compositions. Any information provided by this device is in no way meant to treat, cure or prevent any disease or illness from happening. If in doubt, contact your physician.

This device is contraindicated for any female subject who may be suspected of, or is pregnant. Otherwise, the effects of this device on the fetus are unknown.

This device is contraindicated for any person who is connected to a wearable or implantable electronic device or instrument such as a pacemaker or defibrillator.

This device should not be used for anyone who is acutely or chronically ill because of suffering from a disease or taking medications that affect your water levels. The accuracy of readings for these patients has not been verified. Specific medical advice should be obtained from a physician.

Manufacturer will make available on request circuit diagrams, component parts list, etc.

WARNING: No modifications of this equipment is allowed. This may result in increased EMISSIONS or decreased IMMUNITY of BF-939.

Please use the device according to the user manual. Any misuse can cause electric shock, burns, fire and other unexpected hazards.

Please use and storage the device under the environment which was provided in the user manual.

Don't expose the device to extreme temperatures, direct sunlight, moist or corrosive environment.

Do not step on the scale when your body or feet are wet, especially after bathing or showering to prevent slipping.

INDICATIONS FOR USE

- CAS Body Fat Analyzer measures weight and uses bio-electrical impedance analysis(BIA) technology to estimate body fat, total body water percentage, bone mass, and muscle mass in generally over 18 years old healthy adults.
- It is intended for use in the domestic setting only.


TIPS

To ensure the accuracy of measurement, please follow below instructions when you start measurement.

- Place the scale on a flat, hard surface. Soft surface such as carpet will affect the performance of the scale.
- Step onto the platform with bare feet. Stand still and keep full contact with the electrodes until the measurement is complete.
- Start measurement at least two hours after Getting up or Dinning.
- Avoid measurement immediately after strenuous exercise, sauna or bath, drinking, and dinning.
- Always start measurement in the same time slot and on the same scale located on the same flat, hard surface.
- For maximum accuracy and repeatability , it is recommended that you should use the device in the same time of the day and on the same location.
- The condition of the skin on the bottom of your feet can affect the reading. The natural effects of aging or activity can make this skin hard. Take the reading with clean, slightly damp feet for best accuracy. If you are having a problem on operating this scale please contact customer service.
- Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over-hydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.

DEVICE COMPONENTS


FUNCTION INDICATION

- ◀: Press “◀” to decrease the value in setting mode.
- SET: Touch “SET” to confirm the value in setting mode.
- ▶: Press “▶” to increase the value in setting mode.
- UNIT: Press “UNIT” to turn on the scale.
Press “UNIT” to quit BIA measurement.

LIST

1. Body Fat Analyzer GBF-1257
2. Four AAA-size Batteries (1.5V each)
3. User Manual

LCD DISPLAY

	Body Fat Analysis Result	kg	Kilogram
	Total Body Water Analysis Result	%	Percentage
	Muscle Mass Analysis Result		Man
	Bone Mass Analysis Result		Sportsman
Age	Age		Woman
P8	User ID (Range from P1 to P8)		Sportswoman
cm	Centimeter		Low Battery

GENERAL INSTRUCTIONS

CAS Body Fat Analyzer GBF-1259 applies BIA (Bio-impedance Analysis) technology. A small amount of weak current flows through the human body so as to detect the bio-impedance and estimate body fat, body water, muscle mass and bone mass. The electrical current is small and may not be felt.


This BIA technology is cheap, safe, non-invasive, toxic-free and harmless. It also possesses the characteristics of simple operation and abundant information.

The current mentioned above is less than 0.5mA. However, please be aware that anyone with an wearable or implantable medical electronic instrument, such as a pacemaker, must avoid using this device.

The intended use of this device is for over 18 years old healthy adults.

INSERT THE BATTERIES

- Open the battery door in the back of the scale.
- Insert the batteries(4 x 1.5V AAA) into the battery compartment according to the polarity indications marked inside the compartment.
 - * The digits "8888" will be shown on the LCD.

- Close the battery door and wait until the digits "0.0kg" are shown on the LCD.

⚠ CAUTION


- When the symbol "🔋" appears, the device will power off in about 4s. Then you shall replace with a new set of batteries. Please replace all four batteries at the same time. Do NOT mix the old batteries with the new one.
- Worn batteries are hazardous waste. Do NOT dispose of them together with the household garbage.
- Please refer to the local ordinances and recycling instructions regarding disposal of the worn batteries and scrapped device.
- If you do not intend to use this unit for a prolonged period of time, it is advisable to remove the batteries before storing.

SET UP YOUR PROFILE

The body fat analyzer GBF-1257 supports multiple users (Up to Four). You may follow below instructions to assign User ID and set up your own profile, including Gender, Stature, and Age.

1. Assigning User ID

- With batteries correctly installed, press "SET" key to enter setting.
- The system will request User ID selection first. As pictured to the right, "P1" blinks. The operator may press the function key ◀ or ▶ to select User ID among P1 to P8.
- Touch "SET" key to confirm User ID.


2. Setting Gender

- After confirming User ID, the system will divert to Gender setting automatically.
- As pictured below, the portrait  blinks. The operator may touch the function key ◀ or ▶ to select Gender (Male/Female/Male Athlete/ Female Athlete).
- Touch "SET" key to confirm Gender.
- NOTE:
Athlete Mode: An athlete is defined as an adult who is involved in intense physical activity of approximate 12 hours per week and who has a resting heart rate of approximately 60 beats per minute or less.


3. Setting Stature

- After confirming Gender, the system will divert to Stature setting automatically.
- As pictured below, the digits "160" blinks. The operator may touch the function key ◀ or ▶ to increase or decrease the numeral. (Range: 100 cm to 200 cm)
- You may touch and hold the function key ◀ or ▶ for fast changing the numeral.
- Touch "SET" key to confirm Stature.

4. Setting Age

- After confirming Stature, the system will divert to Age setting automatically.
- As pictured below, the digits "30" blinks. The operator may touch the function key ◀ or ▶ to increase or decrease the numeral.
(Age Setting Range for Normal Mode: 7 to 99 years old,
Age Setting Range for Athlete Mode: 15 to 99 years old)
- You may touch and hold the function key ◀ or ▶ for fast changing the numeral.
- Touch "SET" key to confirm Age.

5. After confirming the Age, the LCD will display "0.0", then you can start measuring.

6. Repeat procedure for a second user, or to change user details

Note: To update or overwrite the memorised data, follow the same procedure, making changes as required.

INITIALISING YOUR SCALE

1. Press the platform centre and remove your foot.
2. "0.0" will be displayed.
3. The scale will switch off and is now ready for use.
This initialisation process must be repeated if the scale is moved.
At all other times step straight on the scale.


WEIGHT ONLY OPERATION

Your CAS Glass Body Fat Analyzer will operate as a conventional weight - reading scale. No special programming steps are required.

Once the scale is initialized, as previously described, you may simply step on the scale to measure your current weight. For only weight reading,


1. Position the scale on a flat, hard surface. Carpeted or uneven floors may affect accuracy.
2. Step onto the scale platform and remain still while the scale computes your weight.
3. The scale will display your weight value.
4. The scale will automatically turn off after a few seconds.


FIRST MEASUREMENT

- Please process the user setting before your first measurement. (Refer to [Set Up Your Profile](#) for more details)

STEP 1: Step on the platform barefooted.


STEP 2: Stand still and keep full contact with the electrodes until the progress bar stops flashing and the measuring results displays.


STEP 3: Then the measuring results will be displayed sequentially three times: Body Fat, Total Body Water, Muscle Mass, Bone Mass.


- If it fails to complete the analysis, the LCD will only display the weight data. (To find out the solutions, please refer to [Troubleshooting](#) for more details.)

DAILY MEASUREMENT


- With original SENSE ON patent technology, GBF-1257 will automatically switch on as you step on the platform barefooted.
- Stand still and keep full contact with the electrodes until the progress bar stops flashing.

- According to the analysis results, the system will automatically identify the possible User ID with most similar history records. Then the measuring results will be displayed sequentially three times.

- When the system finds out two or more users with similar history records, it will notify you to choose between, for example, P1 and P2. You may press ◀ key for P1 and or ▶ key for P2. The measuring results will then be sorted into the User ID you selected and displayed sequentially three times.


- If it fails to identify the possible User ID, the LCD will only display the weight data. (To find out the solutions, please refer to [Troubleshooting](#) for more details.)


BIA MEASUREMENT


- Press the platform centre and remove your foot or press UNIT button to turn on the scale. "0.0" will be displayed.

- Press the DOWN key shortly three times and then press and hold the DOWN key until the indicator of Body Fat lights up.


- Press SET button to set up your own information, including User ID, Gender, Stature, and Age. (Refer to *Set Up Your Profile* for more details)
- Step on the platform barefooted when the LCD displays "0.0".
- Stand still and keep full contact with the electrodes until the progress bar stops flashing.

- Then the value of BIA will be shown on LCD and then turn off.


- NOTE:
If you press the DOWN key shortly three times and then press and hold the DOWN key, it will quit BIA measurement, and display the measurement results three times: Body Fat, Total Body Water, Muscle Mass, Bone Mass.

BODY FAT INDICATOR

The Progress Bar at the top of the LCD Display also functions as the indicator of body fat level. When the measuring result is displayed after analysis, the varying length of the bar indicates different levels of body fat.

(Please refer to *The body fat ratio* for more details.)

ERROR PROMPT

Error	Description	Solution
	Overload. The device will power off.	Stop using this scale for measurement.
	Low Battery. The device will power off in four seconds.	Replace all four batteries at the same time. Please purchase the authorized batteries for replacement.

WHEN MEASURING ...

Problem	Root Cause	Solution
Abnormal measuring results: - Too high; OR - Too low; OR - Huge difference between two recent measurement.	Incorrect posture	Please step on the platform barefooted and stand still.
	The device is located on the soft ground such as a carpet OR on a rugged surface.	Please place the device on a flat, hard surface.
	Cold body that may results in bad blood circulation.	Warm up your hands and feet to resume blood circulation and then measure again.
	Cold Electrodes.	Place the device in a warm room for a while and then measure again.
	Either your hands or your feet are too dry.	Wipe your feet with a damp cloth, keeping them slightly damp when starting measurement.
No display on LCD when the device powers on.	Batteries not yet installed.	Install the batteries. (Please refer to <i><u>Insert the Batteries</u></i>)
	Worn batteries.	Replace all four batteries at the same time. Please purchase the authorized batteries for replacement.
CANNOT proceed to analyze body fat, total body water, muscle mass.	Step onto the platform wearing socks or shoes.	Please keep barefooted during the measurement, and keep full contact with the electrodes as well.
	The system cannot identify the possible User ID with most similar data.	Please assign a User ID following the instruction in <i><u>Set Up Your Profile</u></i> .
	The user fails to select the User ID from what the system found.	Please assign a User ID following the instruction in <i><u>Set Up Your Profile</u></i> .
The device powers off automatically.	Low battery.	Replace all four batteries at the same time. Please purchase the authorized batteries for replacement.

SPECIFICATIONS

Product Name	Body Fat Analyzer (GBF-1257)
Dimension	Scale: 320x330x27.2mm (Approximately)
Net Weight	Approximately 1.75kg (Excluding the dry cells)
Display	Digital LCD with White Backlight
Measurement Unit	Kilogram
Measurement Range	5kg to 180kg
Division	0.1kg
Accuracy	5-50kg: ± 0.3 kg; 50-100kg: ± 0.4 kg; 100-150kg: ± 0.5 kg; 150-180kg: ± 0.7 kg
Working Environment	Temperature: 5°C ~ 40°C Humidity: $\leq 90\%$ RH
Storage Environment	Temperature: -20°C ~ 60°C Humidity: $\leq 93\%$ RH
Power Source	6V (Four AAA-size Batteries)
Turn on Method	SENSE ON technology
Auto-OFF	The scale will turn off after about 10 seconds if there is no operation.
Accessories	1. Four AAA-size Batteries 2. User Manual
Mode of Operation	Continuous Operation

About the Accuracy of This Product

- This product passes strict inspection before delivery and therefore its accuracy is guaranteed by the manufacturer. Please refer to the above table for the descriptions on accuracy.
- This product is specially designed for body fat analysis as well as weight measurement. It should NOT be used by anyone during the process of transaction for verification of goods' weight.

MAINTENANCE

When carrying out usual maintenance, please ensure practice of the following Do's and Don'ts:

- DO use a dry soft cloth to wipe the dust.
- DO use a wet soft cloth, dipped into water and wrung out, to wipe the dirt. Then use a dry soft cloth to dry up the device.
- DON'T wash the device with water or immerse it in water.
- DON'T use propellant, abrasive or other chemicals to wipe the dirt in avoidance of discolor or malfunction.
- DON'T disassemble this device. If you have any problems, please contact CAS.
(Please refer to *Warranty* for contact information)
- Do not dispose of batteries in fire. Batteries may explode or leak. Remove the batteries if the scale will not be used for a long period.

WARRANTY

- CAS warrants its products free of defects in materials and workmanship in normal use for a period of ONE years from the date of retail purchase.
- This warranty does NOT cover damages caused by misuse or abuse, including but not limited to: Failure caused by unauthorized repairs or modifications; Damage caused by shock or drop during transportation; Failure caused by improper operation inconsistent with the instructions stated in this user manual; Malfunction or damage from failure to provide the recommended maintenance; Damage caused by improper use of power supply.
- Should this device require maintenance (or replacement at our option) under warranty, please deliver the original package to GUANGDONG TRANSTEK MEDICAL ELECTRONICS CO., LTD prepaid. Please return the store receipt (with the retail purchase date) and a note with reasons to return on it as well.

THE BODY FAT RATIO (UNIT: %) : 5%-75%

Age	Female				Female			
	Underfat	Healthy	Overfat	Obese	Underfat	Healthy	Overfat	Obese
<39	<21	21.1-33	33.1-39	>39	<8	8.1-20	20.1-25	>25
40-59	<23	23.1-34	34.1-40	>40	<11	11.1-22	22.1-28	>28
>59	<24	24.1-36	36.1-42	>42	<13	13.1-25	25.1-30	>30

THE BODY MUSCLE RATIO (%): 35%-90%

Age	Female	Male
	Healthy	Healthy
18-99	>34	>40

THE BONE RANGE (KG) ≤ 10KG

Gender	Female			Male		
Weight	<45kg	45-60kg	>60kg	<60kg	60-75kg	75kg
Bone	1.8kg	2.2kg	2.5kg	2.5kg	2.9kg	>3.2kg

HEALTH TIPS - ABOUT BODY FAT

Fat is essential for human body. It can not only store energy and protect viscera, but also regulate body temperature and maintain normal physiological function of human body. However, too much body fat is harmful to human body. It is always accompanied by Fatty Liver, diabetes, coronary heart disease, etc.

Therefore self-measuring and self-monitoring body fat level are beneficial to your health. Since we can't judge body fat level simply by our weight, this body fat analyzer GBF-1257, with BIA (Bio-impedance Analysis) technology applied, is an accurate device that offers a quick and comfortable way to obtain your body fat level.

EMC GUIDANCE

Table 1 – Guidance and MANUFACTURER'S declaration
– ELECTROMAGNETIC EMISSIONS for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 2	The device must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class B	
Harmonic missions IEC 61000-3-2	Not applicable	
Voltage fuctuations /flicker emissions IEC 61000-3-3	Not applicable	

Table 2 – Guidance and MANUFACTURER'S declaration – electromagnetic IMMUNITY
 – for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

**Table 4 – Guidance and MANUFACTURER’S declaration – electromagnetic IMMUNITY
– for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

Guidance and manufacturer’s declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>Not applicable</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d = 3,5/\sqrt{V_1} \sqrt{P}$ $d = 1,167 \sqrt{P}$ 80MHz to 800 MHz $d = 2,333 \sqrt{P}$ 800MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.</p>			
<p>b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than $[V_1]$ V/m.</p>			

Table 6 – Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM
 – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the device			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter		
	150 kHz - 80 MHz $d = 3,5/\sqrt{P}$	80 MHz - 800 MHz $d = 1,167 \sqrt{P}$	800MHz to 2,5 GHz $d = 2,333 \sqrt{P}$
0,01	Not applicable	0.117	0.233
0,1	Not applicable	0.369	0.738
1	Not applicable	1.167	2.333
10	Not applicable	3.690	2.333
100	Not applicable	11.67	23.33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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