KD-595

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

Cuff tube plug

Rubber tube

Blood pressure level
classification indicator

Memory

Start

LCD

Function Description

Based on oscillometric method and silicon integrate pressure sensor technology, blood pressure and heart rate can be measured automatically and non-invasively with voice function which helps to ease the operation. New blood pressure classification indicator visually indicates user's blood pressure category, 60 memories with date and time display. Big LCD screen will show the result of measurement and irregular heartbeat (if any) .



Specification

- | | |
|---|--|
| (1) Model: KD-595 | (9) Pressure measuring range:0-300mmHg |
| (2) Classification: Class II, Type B | (10) Measuring accuracy: ± 3 mmHg |
| (3) Machine size: 100mmx120mmx45mm | (11) Heart rate range: 30-180/min |
| (4) Cuff size: 520mmx140mm | (12) Environmental temperature for using: 5°C~40°C |
| (5) Weight: about 300g (exclude batteries) | (13) Environmental humidity for using: <85% |
| (6) Measuring method: oscillometric method, automatic air inflation and measurement | (14) Environmental temperature for storage: -20°C~55°C |
| (7) Memory volume: 60 times with time and date displaying | (15) Environmental humidity for storage: <95% |
| (8) Power source: DC 6V or AAx4 batteries | (16) Environmental pressure: atmosphere pressure |
| | (17) Battery life: Approx.2 months with 3-min. usage per day |

Measuring Procedures

1. Battery Loading

- a. Open battery cover at the back of the machine.
 - b. Load four "AA" size batteries.
Please pay attention to polarity.
 - c. Close the battery cover.
- ★ When LCD shows battery signal it means you must change the batteries at once. Rechargeable batteries with a voltage of 1.2V are not suitable for this product.
 - ★ Take out all the batteries if the device will be not used for a long time to avoid battery leakage and relevant damage.


External power supply socket

The equipment, the batteries and cuff used in it can be disposed as industrial rubbish at the end of their usage or according to local regulations.

2. Clock Adjusting


- a. When the machine is turned off, press the button "START" and "MEM" at the same time, the number of month will twinkle on LCD. Press button "START" continuously, the numbers of month, date, hour, and minute will twinkle in sequence. While the number is twinkling, press button "MEM" and the corresponding number will increase by one. Hold on pressing the button "MEM" the number will increase fast.
- b. After pressing "START" button the fifth time, numbers on LCD will not twinkle, then you exit the time adjusting mode. See picture1.
- c. After you change the batteries, you have to readjust the date and time.

picture 1


3. Voice Setting

- a. Voice language setting: In clock display mode, you can select the voice language by keeping on press the button "MEM". Now LCD twinkles "L0", "L1", "L2"... "Ln"(n means the number of language) circularly. "L0" represents closing voice function, "L1" represents language 1, "L2" represents language 2, ..."Ln" represents language n and so on. The monitor will speak out the corresponding language name of the displayed language code. You can select the wanted language by releasing button "MEM" when display the corresponding language code. See picture 2.
- b. If the number of language is only one, LCD will twinkle "ON" and "--" circularly. Here "ON" represents opening the voice function and "--" represents closing the voice function.
- c. Voice volume setting: You can select the voice volume by keeping on press the button "MEM" again after language setting. Now the bars on the LCD

picture 2



change and the number of the bars indicates the volume. You can select the wanted volume by releasing button "MEM" when display the corresponding bars. See picture 3.

picture 3

4. Cuff Connecting

- Be sure there is no leakage when placing the rubber tube plug of cuff socket.
- Put the end of cuff through metal hoop. Make sure that the two nylon outlets are located on the outside surface of the cuff.

5. Cuff Swathing

- Place cuff 1-2cm above elbow joint.
 - ★ Place palm upside and let it be parallel with rubber tube.
 - ★ Measuring at the same arm each time.
 - ★ Measuring on the bare arm, especially for the feeble pulse.
 - ★ When the cuff is dirty, detach it from the equipment, wash the cuff by hand with proper detergent and rinse it in the enough cold water. Dry in air. Never iron it.

- b. Turn the cuff end outward and tighten it by pulling through the metal loop. Then press the loop onto the surface of cuff.
- c. Do not place the cuff too loose or too tight; spacing of a finger is appreciated.


Sitting Measurement

- a. Sit upright
- b. Place arm horizontally on the chair or table with the palm upside.
- c. Cuff position should be at the same level as the heart.

Lying down Measurement

- a. Lie on your back
- b. Naturally place arm and hand straight, with palm upside.
- c. Keep the cuff at the same level as the heart.

- ★ During measurement, please make sure the cuff stays at the same height as the heart. Either above or below the position, the accuracy will be affected.
- ★ Stay calm for 5 minutes before measuring and keep yourself comfortable, steady and quiet during measuring.
- ★ Do not sway your body or move the monitor or press the rubber tube. Avoid swaying your arms.

6. Measuring Process

- a. After your body is in a comfortable position, press the "START" button. After a long "beep", BPM will verify itself, and LCD will twinkle automatically. See picture 4.


picture 4

- b. If it already has memory, then LCD will show the latest result. See picture 5. If machine does not have memory before, the LCD will show five "0", while it verify itself. See picture 6.


picture 5

picture 6

- c. If the voice function is on, the monitor will speak out the measurement tips with that language.

- d. After two seconds (maybe longer because of verifying itself), LCD shows one "0".
- e. If the voice function is on, the monitor will speak out the prompt that the measurement will begin.
- f. Then the machine begin to pressurize automatically.
- g. After measurement, the result of systolic, diastolic and pulse rate will be all displayed on LCD. And irregular heartbeat symbol "⚡❤️" (if any) and blood pressure classification indicator "■" will twinkle, see picture 7. The result will be automatically memorized by itself.
- h. If the voice function is on, the monitor will speak out the displayed blood pressure, heart rate and blood pressure level according to the classification standard of WHO(World Health Organization).
- i. Pressing "START" button after measuring, machine will be turned off. It will be automatically power off after 1 minute of none use.


picture 7

7. Assessing High Blood Pressure for Adults

The follow standards for assessing high blood pressure (without regard to age or gender) have been established as a guide line according to WHO (World Health Organization) standard. Please note that other risk factors (e.g. diabetes, obesity, smoking, etc.) need to be taken into consideration and may affect these figures. Consult with your physician for accurate assessment.

Classification of blood pressure for adults

BLOOD PRESSURE CLASSIFICATION	SBP mmHg	DBP mmHg	COLOR INDICATOR
Optimal	<120	<80	GREEN
Normal	120-129	80-84	GREEN
High-Normal	130-139	85-89	GREEN
Stage 1 Hypertension	140-159	90-99	YELLOW
Stage 2 Hypertension	160-179	100-109	ORANGE
Stage 3 Hypertension	≥ 180	≥ 110	RED

WHO/ISH Definitions and Classification of Blood Pressure Levels

8. Reading Result of Memories

- a. Under the clock-displaying mode, press button "MEM" and enter the mode of reading results. LCD will show how many results it has, then the latest result will be displayed, which is 01 for memory number, see picture 8. Continuously press the "MEM" button, the memory number will increase by one, see picture 9.
- b. If the voice function is on, the monitor will speak out the displayed memory sequence number, blood pressure, heart rate and blood pressure level according to the classification standard of WHO(World Health Organization).


picture 8


picture 9

- c. Press "MEM" button again after the last memory, the machine will power off. Under the memory-displaying mode, it will be automatically power off after ten seconds of none use. Or you can press "START" button to shut down the machine.
- d. If it has no memory, then LCD will show "0" for systolic, diastolic, and pulse rate, and show "00" for the number of memory. While press the "MEM" or "START", machine will power off.

9. Deletion of Memories

Under the mode of reading memories, holding on pressing "MEM" button for three seconds, all memories will be deleted after three "beep". LCD will show the picture 10 for one second, then come to the mode without memory, see picture 11.


picture 10

picture 11

10. Shut Down

After measurement, press button "START" to turn off the machine; the machine will be automatically power off after 1 minute of none use.

Abnormality & Troubleshooting (1)

Abnormality	Reason	Checkout
LCD shows battery symbol 	Batteries low	Change new batteries
Shows abnormal result	Cuff not tightened properly or its position is incorrect	Tighten cuff correctly and refer to "Cuff Swathing"
	The arm is moved during measuring or BPM is moved	Stay calm, remains steady. Do not move BPM
	Irregular heartbeat	Test again for light irregular heartbeat patients. It is inappropriate for serious irregular heartbeat patients to use this machine.

Abnormality & Troubleshooting (2)

Abnormality	Reason	Checkout
Shows abnormal result	Speaking, frightened nervous or excited during measurement	Do not speak, take deep breath 2~3 times to relax yourself
	Wrong position	Adjust position; refer to "Cuff Swathing"
	Some interference in inflation or wrong operation during measuring	Refer to the inflation step in "Measuring Process"
LCD shows "EE"	Insufficient inflation	Measure again
	Some interference in inflation or wrong operation during measuring	Measure again
BPM no response when you press button, BPM hang	Incorrect operation	Take out batteries for a while, and then reinstall batteries.

Note: If you cannot resolve the problem, you can contact factory or its designated service department.

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