HEALTHCARE SOLUTION™

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Thank you very much for selecting the CAS Body Fat Analyzer GBF-1039.

Please do read the user manual carefully and thoroughly so as to ensure the safe usage of this product, and keep the manual well for further reference in case you have problems.

TABLE OF CONTENT

Safety Information	
- Safety and Usage Information	24
- Tips	25
Overview	
- Device Components	26
- List	26
- LCD Display	27
Initial Start-Up	
- General Instructions	28
- Insert the Batteries	28
- Set Up Your Profile	29
- Initialising Your Scale	31
Start Measuring	
- Weight Only Operation	32
- First Measurement	33
- Daily Measurement	34
Manage Your Health	
- Data Enquiry In Your Scale	35
Troubleshooting	
- Error Prompt	36
- When Measuring	36
Specifications	38
Maintenance	39
Warranty	39
Appendix	
- Health Tips-About Body Fat	40
- Health Tips-About Body Water	41
- Health Tips-About Muscle Mass	42
- EMC Guidance	42

Safety Information

Safety and Usage Information

The warning signs and symbols are essential to ensure your correct and safe use of this product and protect you and others from injury. Please kindly find the meanings of the warning signs and symbols, which you may encounter in the label and user manual, as follows:

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "MANUFACTURER"
	Symbol for "DIRECT CURRENT"		Symbol for "COMPLIES WITH EU REQUIREMENTS"
	Symbol for "MANUFACTURE DATE"		Symbol for "ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
SN	Symbol for "SERIAL NUMBER"		

CAUTION

- CAS's Body Fat Analyzer GBF-1039 offers you a seamless way to manage your health. Please be aware that this device is designed for healthy population over 10 years old self-measuring and self-monitoring body compositions. Any information provided by this device is in no way meant to treat, cure or prevent any disease or illness from happening. If in doubt, contact your physician.
- This device is contraindicated for any female subject who may be suspected of, or is pregnant. Otherwise, the effects of this device on the fetus are unknown.
- This device is contraindicated for any person who is connected to a wearable or implantable electronic device or instrument such as a pacemaker or defibrillator.
- This device should not be used for anyone who is acutely or chronically ill because of suffering from a disease or taking medications that affect your water levels. The accuracy of readings for these patients has not been verified. Specific medical advice should be obtained from a physician.
- Manufacturer will make available on request circuit diagrams, component parts list, etc.

WARNING:

No modifications of this equipment is allowed. This may result in increased EMISSIONS or decreased IMMUNITY of GBF-1039. Please use the device according to the user manual. Any misuse can cause electric shock, burns, fire and other unexpected hazards. Please use and storage the device under the environment which was provided in the user manual. Don't expose the device to extreme temperatures, direct sunlight, moist or corrosive environment. Don't step on the scale when your body or feet are wet, especially after bathing or showering to prevent slipping. Please keep the device out of reach of infants, children or pets, since inhalation or swallowing of small parts is dangerous or even fatal. Please don't use the device out of lifetime any more because degraded sensors and electrodes, or loosened electrodes can degrade performance or cause other problems. If you have an allergy to stainless steel, please avoid to contact to the electrodes of the device. The patient is an intended operator. The patient can measure, transmit data and charge battery under normal circumstances and maintain the device and its accessories according to the user manual.

INDICATIONS FOR USE

- The CAS Body Fat Analyzer measures weight and uses bio-electrical impedance analysis (BIA) technology to estimate body fat, total body water percentage, and muscle mass in generally healthy children 10-17 years old and healthy adults.
- It is intended for use in the domestic setting only.



Tips

To ensure the accuracy of measurement, please follow below instructions when you start measurement.

- Place the scale on a flat, hard surface. Soft surface such as carpet will affect the performance of the scale.
- Step onto the platform with bare feet. Stand still and keep full contact with the electrodes until the measurement is complete.
- Start measurement at least two hours after Getting up or Dinning.
- Avoid measurement immediately after strenuous exercise, sauna or bath, drinking, and dinning.
- Always start measurement in the same time slot and on the same scale located on the same flat, hard surface.
- For maximum accuracy and repeatability , it is recommended that you should use the device in the same time of the day and on the same location.
- The condition of the skin on the bottom of your feet can affect the reading. The natural effects of aging or activity can make this skin hard. Take the reading with clean, slightly damp feet for best accuracy. If you are having a problem on operating this scale please contact customer service.
- Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over-hydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.

Overview

Device Components

List

1. Body Fat Analyzer GBF-1039
2. Two CR2032 Batteries (3V each)
3. User Manual

LCD Display

fat	Body Fat Analysis Result	cm	Centimeter
TBW	Total Body Water Analysis Result	kg	Kilogram
MUS	Muscle Mass Analysis Result	%	Percentage
Underfat	Body Fat Level: Underfat	age	Age
Healthy	Body Fat Level: Healthy		Male
Overfat	Body Fat Level: Overfat		Female Athlete
Obese	Body Fat Level: Obese		Female
8	User ID (Range from 1 to 4)		Female Athlete

Initial Start-Up

General Instructions

CAS Body Fat Analyzer GBF-1039 applies BIA (Bio-impedance Analysis) technology. A small amount of weak current flows through the human body so as to detect the bio-impedance and estimate body fat, body water, muscle mass. The electrical current is small and may not be felt.

This BIA technology is cheap, safe, non-invasive, toxic-free and harmless.

It also possesses the characteristics of simple operation and abundant information.

The current mentioned above is less than 0.5mA. However, please be aware that anyone with an wearable or implantable medical electronic instrument, such as a pacemaker, must avoid using this device.

The intended use of this device is for healthy children 10-17 years old and healthy adults.

Insert the Batteries

- Open the battery door in the back of the scale.
- Insert the batteries (2 x CR2032) into the battery compartment according to the polarity indications marked inside the compartment.
- Close the battery door and wait until the digits "0.0kg" are shown on the LCD.

CAUTION

- When the symbol "Lo" appears, the device will power off in about 4s. Then you shall replace with a new set of batteries. Please replace both two batteries at the same time. Do NOT mix the old batteries with the new one.
- Worn batteries are hazardous waste. Do NOT dispose of them together with the household garbage.
- Please refer to the local ordinances and recycling instructions regarding disposal of the worn batteries and scrapped device.
- If you do not intend to use this unit for a prolonged period of time, it is advisable to remove the batteries before storing.

Set Up Your Profile

The body fat analyzer GBF-1039 supports multiple users (Up to Four). You may follow below instructions to assign User ID and set up your own profile, including Gender, Stature, and Age.

- Turn on Bluetooth and the App. Make sure both are ON when pair-up is proceeding.
- Press and hold "UNIT" button in the back of the scale to start pair-up.

1. Assigning User ID

- With batteries correctly installed, touch "SET" key to enter setting.
- The system will request User ID selection first. As pictured below, "1" blinks.
The operator may touch the function key < or > to select User ID among 1 to 4.
- Touch "SET" key to confirm User ID.

2. Setting Gender

- After confirming User ID, the system will divert to Gender setting automatically.
As pictured below, the portrait  blinks. The operator may touch the function key < or > to select Gender (Male/Female/Male Athlete/ Female Athlete).

3. Setting Stature

- After confirming Gender, the system will divert to Stature setting automatically.
As pictured below, the digits "165" blinks. The operator may touch the function key < or > to increase or decrease the numeral.
- You may touch and hold the function key < or > for fast changing the numeral.
- Touch "SET" key to confirm Stature.

4. Setting Age

- After confirming Stature, the system will divert to Age setting automatically.
As pictured below, the digits "30" blinks. The operator may touch the function key < or > to increase or decrease the numeral.
- You may touch and hold the function key < or > for fast changing the numeral.
- Touch "SET" key to confirm Age.

5. After confirming the Age, the LCD will display "0.0", then you can start measuring.

6. Repeat procedure for a second user, or to change user details.

Note: To update or overwrite the memorised data, follow the same procedure, making changes as required.

Initialising Your Scale

1. Press the platform centre and remove your foot.
2. "0.0" will be displayed.

3. The scale will switch off and is now ready for use.
This initialisation process must be repeated if the scale is moved.
At all other times step straight on the scale.


Start Measuring

Weight Only Operation

Your CAS Body Fat Analyzer will operate as a conventional weight - reading scale. No special programming steps are required. Once the scale is initialized, as previously described, you may simply step on the scale to measure your current weight.

For only weight reading.

1. Position the scale on a flat, hard surface. Carpeted or uneven floors may affect accuracy.

2. Step onto the scale platform and remain still until the LCD stops displaying moving "o".

3. The User ID blinks (eg, 1), you can select User ID and set up your own profile, and then start measuring when LCD displays 0.0.
(Refer to Set Up Your Profile for more details)


4. The scale will automatically turn off after a few seconds if there is no operation.

First Measurement

- Please process the user setting before your first measurement.
(Refer to **Set Up Your Profile** for more details)
- Position scale on a firm flat surface.
- Touch SET button.
- While the user number of the latest measurement is flashing, select your user number by touching < or > button. After 3s, the selected number is locked, the scale shows zero reading.
(Note: When the user number is flashing, if you touch SET button again, it will enter the setting mode, after all the settings are finished, it will display zero reading.)

STEP 1 : Step on the platform barefooted


STEP 2 : Stand still and keep full contact with the electrodes until the LCD stops displaying moving "o".

STEP 3 : Then the measuring results will be displayed sequentially three times:

Body Fat, Total Body Water, Muscle Mass, Weight.

And the personal information will be displayed once, and then turn off.


* If it fails to complete the analysis, the LCD will only display the weight data.
(To find out the solutions, please refer to **Troubleshooting** for more details.)

Daily Measurement

With original SENSE ON patent technology, GBF-1039 will automatically switch on as you step on the platform barefooted.

According to the analysis results, the system will automatically identify the possible User ID with most similar history records. Then the measuring results will be displayed sequentially three times.

If the system fails to identify the possible User ID, the digits 1234 will be shown on LCD, then you may touch < or > to select the User ID manually. The measuring results will then be sorted into the User ID you selected and displayed sequentially three times.

If it fails to identify the possible User ID, the LCD will only display the weight data. (To find out the solutions, please refer to **Troubleshooting** for more details.)

Manage Your Health

Data Enquiry in Your Scale

When the scale display "0.0", touch the function key < or > to check the latest analysis record stored in the device.

It will start from the User ID that completes the last measurement.

You may touch the function key < or > to view other users' latest analysis records.

The analysis result will be displayed once in the following order:

Body Fat, Total Body Water, Muscle Mass.

When there are no pieces of measuring result stored under a specific User ID, for example, 3, the LCD will display " - - - - " as pictured below.

After displaying the analysis result, the device will divert to the weighing interface automatically (as pictured below). You may start measuring and the measuring result will be recorded under the referring User ID.

Troubleshooting

Error Prompt

Error	Description	Solution
	Overload. The device will power off.	Stop using this scale for measurement.
	Low Battery. The device will power off in four seconds.	Replace both two batteries at the same time. Please purchase the authorized batteries for replacement.
	Measuring error.	Step off the platform, and measure again.

When Measuring

Problem	Root Cause	Solution
Abnormal measuring results: - Too high; OR - Too low; OR - Huge difference between two recent measurement.	Incorrect posture	Please step on the platform barefooted and stand still.
	The device is located on the soft ground such as a carpet OR on a rugged surface.	Please place the device on a flat, hard surface.
	Cold body that may results in bad blood circulation.	Warm up your hands and feet to resume blood circulation and then measure again.

Problem	Root Cause	Solution
Abnormal measuring results: - Too high; OR - Too low; OR - Huge difference between two recent measurement.	Cold Electrodes.	Place the device in a warm room for a while and then measure again.
	Either your hands or your feet are too dry.	Wipe your feet with a damp cloth, keeping them slightly damp when starting measurement.
No display on LCD when the device powers on.	Batteries not yet installed	Install the batteries. (Please refer to Insert the Batteries)
	Worn batteries	Replace both two batteries at the same time. Please purchase the authorized batteries for replacement.
CANNOT proceed to analyze body fat, total body water, muscle mass.	Step onto the platform wearing socks or shoes.	Please keep barefooted during the measurement, and keep full contact with the electrodes as well.
	The system cannot identify the possible User ID with most similar data.	Please assign a User ID following the instruction in Set Up Your Profile.
	The user fails to select the User ID from what the system found	Please assign a User ID following the instruction in Set Up Your Profile.
The device powers off automatically.	Low battery.	Replace both two batteries at the same time. Please purchase the authorized batteries for replacement.

Specifications

Product Name	Body Fat Analyzer (GBF-1039)
Dimension	Scale: 298.9x296.6x32.4mm(Approximately)
Net Weight	Approximately 1.48kg (Excluding the cells)
Display	Digital LCD
Measurement Unit	Kilogram
Measurement Range	5kg to 150kg
Division	0.1kg
Accuracy	5-50kg: ± 0.3 kg; 50-100kg: ± 0.4 kg; 100-150kg: ± 0.5 kg;
Working Environment	Temperature: 5°C to 40°C Humidity: 15% RH to 93%RH Atmospheric pressure: 70kPa to 106kPa
Storage Environment	Temperature: -20°C to 60°C Humidity: $\leq 93\%$ RH Atmospheric pressure: 50kPa to 106kPa
Power Source	6V (Two CR2032 Batteries)
Turn on Method	SENSE ON technology
Auto-OFF	The scale will turn off after about 10 seconds if there is no operation.
Accessories	1. Two CR2032 Batteries 2. User Manual
Mode of Operation	Continuous Operation

About the Accuracy of This Product

- This product passes strict inspection before delivery and therefore its accuracy is guaranteed by the manufacturer.
Please refer to the above table for the descriptions on accuracy.
- This product is specially designed for body fat analysis as well as weight measurement. It should NOT be used by anyone during the process of transaction for verification of goods' weight.

Maintenance

Maintenance

When carrying out usual maintenance, please ensure practice of the following Do's and Don'ts :

- DO use a dry soft cloth to wipe the dust.
- DO use a wet soft cloth, dipped into water and wrung out, to wipe the dirt. Then use a dry soft cloth to dry up the device.
- DON'T wash the device with water or immerse it in water.
- DON'T use propellant, abrasive or other chemicals to wipe the dirt in avoidance of discolor or malfunction.
- DON'T disassemble this device. If you have any problems, please contact CAS. (Please refer to Warranty for contact information)
- Do not dispose of batteries in fire. Batteries may explode or leak. Remove the batteries if the scale will not be used for a long period.

Warranty

- CAS warrants its products free of defects in materials and workmanship in normal use for a period of ONE year from the date of retail purchase.
- This warranty does NOT cover damages caused by misuse or abuse, including but not limited to:
 - Failure caused by unauthorized repairs or modifications;
 - Damage caused by shock or drop during transportation;
 - Failure caused by improper operation inconsistent with the instructions stated in this user manual;
 - Malfunction or damage from failure to provide the recommended maintenance;
 - Damage caused by improper use of power supply.
- Should this device require maintenance (or replacement at our option) under warranty, please deliver the original package to CAS prepaid. Please return the store receipt (with the retail purchase date) and a note with reasons to return on it as well.

Appendix

Health Tips - About Body Fat

Human body is made up of, amongst other things, a percentage of fat. Body fat is vital for a healthy, functioning body, protects vital organs, helps regulate body temperature, stores vitamins and helps the body sustain itself when food is scarce. However, too much body fat or indeed too little body fat will damage to your health. It is difficult to gauge how much body fat we have in our bodies simply by looking at ourselves in the mirror.

This is why it is important to measure and monitor your body fat percentage. Body fat percentage gives you a better measure of fitness than weight alone—the composition of your weight loss could mean you are losing muscle mass rather than fat—you could still have a high percentage of fat even when a scale indicates 'normal weight'.

The table as follows may be used as a guide:

The body fat ratio (Unit: %)

Age	Female				Male			
	Underfat	Healthy	Overfat	Obese	Underfat	Healthy	Overfat	Obese
<39	<21	21.1-33	33.1-39	>39	<8	8.1-20	20.1-25	>25
40-59	<23	23.1-34	34.1-40	>40	<11	11.1-22	22.1-28	>28
>59	<24	24.1-36	36.1-42	>42	<13	13.1-25	25.1-30	>30

Health Tips - About Body Water

Body water is the single most important component of body weight. It represents over half of your total weight and almost two thirds of your lean body mass (predominantly muscle). Water performs a number of important roles in the body:

All the cells in the body, whether in the skin, glands, muscles, brain or anywhere else, can only function properly if they have enough water. Water also plays a vital part in regulating the body's temperature balance, particularly through perspiration.

The combination of your weight and fat measurement could appear to be 'normal' but your body hydration level could be insufficient for healthy living.

The table as follows may be used as a guide:

The body water ratio (Unit: %)

Age	Female			Male		
	Low	Normal	High	Low	Normal	High
10-15	<57	57.1-67	>67	<58	58.1-72	>72
16-30	<47	47.1-57	>57	<53	53.1-67	>67
31-60	<42	42.1-52	>52	<47	47.1-61	>61
61-85	<37	37.1-47	>47	<42	42.1-56	>56

Health Tips - About Muscle Mass

According to the American College of Sports Medicine (ACSM), lean muscle mass may decrease by nearly 50 percent between the age of 20 and 90. If you do nothing with it you're losing muscle and increasing fat. It is also important to know your muscle mass % during weight reduction. At rest, the body burns approximately 110 additional calories for each kilo of muscle gained. Benefits of gaining muscle mass include:

- Reversing the decline in strength, bone density and muscle mass with aging
- Maintenance of flexible joints
- Guide weight reduction when combined with a healthy diet.

The table as follows may be used as a guide:

The body muscle ratio (%):

Age	Female	Male
	Healthy	Healthy
10-85	>34	>40

EMC Guidance

1. MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.
2. Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance $d = 3,3$ m away from the equipment.

(Note. As indicated in Table 6 of IEC 60601-1-2:2007 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2 W yields $d = 3,3$ m at an IMMUNITY LEVEL of 3 V/m)

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